DANSK

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**ARTIQ**

Analyzed: 2022-11-21 15:55 UTC±00:00  
 Report ID: DEMO  
 Age: 73    Gender: Female    Current Smoker: No    Pack-Years: 1

**Interpretation of lung function tests**  
 Moderate obstructive spirometry.  
 No significant reversibility after postbronchodilator test. Change in FEV1 of 90ml or 8.9%, and FVC of 130ml or 6.2%.  
 Mild reduction of diffusion capacity.

**Disease probability**

Disease	Probability
Asthma	17.8%
COPD	56.8%
OBD	22.6%

**Conclusions and suggestions**  
 Highest disease probability based on lung function: COPD

In addition to patient history and clinical exam, the following tests may help determine a final diagnosis:  
 Perform a lung volume test.

**Warnings**  
 No data available on lung volumes. Calculated probabilities are more accurate when full lung function data is provided.

**Legend**  
 COPD: Chronic Obstructive Pulmonary Disease  
 OBD: Other Obstructive Diseases (including: cystic fibrosis, bronchiectasis, bronchidilia)  
 Normal: Normal lung function  
 ILD: Interstitial lung disease (including: idiopathic pulmonary fibrosis, nonspecific interstitial pneumonitis, sarcoidosis)  
 NMD: Neuromuscular disease (including: paralysis of the diaphragm, poliomyelitis, myopathy)  
 PVD: Pulmonary vascular disease (including: pulmonary hypertension, embolism, vasculitis)  
 TD: Thoracic deformity / Pleural disease (including: pneumothorax, lobectomy, chest wall problems, kyphoscoliosis)

This report is approved for clinical use in the EU.  
 Automatically generated by ArtiQ.PFT TEST-5bb4c8b • Manufactured by ArtiQ NV • Leuven, Belgium

info@ArtiQ.eu  
 www.ArtiQ.eu

**1. Interpretation: PFT description as dictated by the international standards**

**2. Analysis: disease probability estimation with machine learning engine**

**3. Decision support: most likely disease**

**4. Further suggestions: best medical practice**

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