DANSK

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Analyzed: 2022-11-21 15:55 UTC±00:00
Report ID: DEMO
Age: 73 **Gender:** Female **Current Smoker:** No **Pack-Years:** 1

<p>Interpretation of lung function tests Moderate obstructive spirometry. No significant reversibility after postbronchodilator test. Change in FEV1 of 90ml or 8.9%, and FVC of 130ml or 6.2%. Mild reduction of diffusion capacity.</p>	<p>1. Interpretation: PFT description as dictated by the international standards</p>			
<p>Disease probability</p> <div style="display: flex; align-items: center;"> <div style="flex: 1;"> <table style="width: 100%; text-align: center;"> <tr> <td style="width: 33%;">Asthma 17.8%</td> <td style="width: 33%;">COPD 56.8%</td> <td style="width: 33%;">OBD 22.6%</td> </tr> </table> </div> <div style="flex: 0.5; font-size: 0.8em; margin-left: 10px;"> <ul style="list-style-type: none"> <input type="checkbox"/> Asthma <input type="checkbox"/> COPD <input type="checkbox"/> OBD <input type="checkbox"/> Normal <input type="checkbox"/> ILD <input type="checkbox"/> NMD <input type="checkbox"/> PVD <input type="checkbox"/> TD </div> </div>	Asthma 17.8%	COPD 56.8%	OBD 22.6%	<p>2. Analysis: disease probability estimation with machine learning engine</p>
Asthma 17.8%	COPD 56.8%	OBD 22.6%		
<p>Conclusions and suggestions Highest disease probability based on lung function: COPD</p>	<p>3. Decision support: most likely disease</p>			
<p>In addition to patient history and clinical exam, the following tests may help determine a final diagnosis: Perform a lung volume test.</p>	<p>4. Further suggestions: best medical practice</p>			
<p>Warnings No data available on lung volumes. Calculated probabilities are more accurate when full lung function data is provided.</p>	<p>5. Warnings</p>			

Legend

COPD Chronic Obstructive Pulmonary Disease
 OBD Other Obstructive Diseases (including cystic fibrosis, bronchiectasis, bronchiolitis)
 Normal Normal lung function
 ILD Interstitial lung disease (including idiopathic pulmonary fibrosis, nonspecific interstitial pneumonitis, sarcoidosis)
 NMD Neuromuscular disease (including paralysis of the diaphragm, poliomyelitis, myopathy)
 PVD Pulmonary vascular disease (including pulmonary hypertension, embolism, vasculitis)
 TD Thoracic deformity / Pleural disease (including: pneumothorax, lobectomy, chest wall problems, kyphoscoliosis)

This report is approved for clinical use in the EU.
 Automatically generated by ArtiQ.PFT TEST-5bb4c8b • Manufactured by ArtiQ NV • Leuven, Belgium

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QUNIQUE GmbH, Bahnhofweg 17, 5610 Wohlen, Schweiz