DANSK

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Analyzed: 2022-11-21 15:55 UTC±00:00
 Report ID: DEMO
 Age: 73 Gender: Female Current Smoker: No Pack-Years: 1

Interpretation of lung function tests
 Moderate obstructive spirometry.
 No significant reversibility after postbronchodilator test. Change in FEV1 of 90ml or 8.9%, and FVC of 130ml or 6.2%.
 Mild reduction of diffusion capacity.

Disease probability

Disease	Probability
Asthma	17.8%
COPD	56.8%
OBID	22.6%

Legend: Asthma, COPD, OBID, Normal, ELD, NMD, PVD, TD

Conclusions and suggestions
 Highest disease probability based on lung function: COPD

In addition to patient history and clinical exam, the following tests may help determine a final diagnosis:
 Perform a lung volume test.

Warnings
 No data available on lung volumes. Calculated probabilities are more accurate when full lung function data is provided.

Legend
 COPD Chronic Obstructive Pulmonary Disease
 OBID Other Obstructive Diseases (including cystic fibrosis, bronchiectasis, bronchiolitis)
 Normal Normal lung function
 ELD Interstitial lung disease (including idiopathic pulmonary fibrosis, sarcoidosis, hypersensitivity pneumonitis, sarcoidosis)
 NMD Neuromuscular disease (including paralysis of the diaphragm, polymyositis, myopathy)
 PVD Pulmonary vascular disease (including pulmonary hypertension, embolism, vasculitis)
 TD Thoracic deformity/ Pleural disease (including pneumectomy, lobectomy, chest wall problems, kyphoscoliosis)

This report is approved for clinical use in the EU.
 Automatically generated by ArtiQ.PFT 1.9.1-36645b • Manufactured by ArtiQ NV • Leuven, Belgium info@ArtiQ.eu www.ArtiQ.eu

1. Interpretation: PFT description as dictated by the international standards

2. Analysis: disease probability estimation with machine learning engine

3. Decision support: most likely disease

4. Further suggestions: best medical practice

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QUNIQUE GmbH, Bahnhofweg 17, 5610 Wohlen, Schweiz