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|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                    |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| <p><b>ARTIQ</b></p> <p>Analyzed: 2025-10-27 13:11<br/>Report ID: 699f20cd-9c3f-4508-a34a-9750adff5d60<br/>Age: 64      Sex at birth: Female      Current Smoker: Yes      Pack-Years: 34</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                                    |
| <p><b>Quality of lung function tests</b><br/>FEV1 quality grade is A. FVC quality grade is B.<br/>The test quality is good and results can be interpreted with confidence.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 1. <b>Quality:</b> Quality of spirometry and meaning                                               |
| <p><b>Physiological interpretation of lung function tests</b><br/>Moderate obstructive lung function.<br/>Bronchodilator response test is not performed.<br/>Hyperinflation and/or air trapping.<br/>Increased airway resistance.<br/>Moderate reduction of diffusion capacity.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | 2. <b>Physiological interpretation:</b> spirometry description as per international standards      |
| <p><b>AI-supported interpretation</b></p> <div> <p>COPD 67.8%</p> <ul style="list-style-type: none"> <li>Asthma</li> <li>COPD</li> <li>OBD</li> <li>Normal</li> <li>ILD</li> <li>NMD</li> <li>PVD</li> <li>TD</li> </ul> </div>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 3. <b>AI-supported interpretation:</b> disease probability estimation with machine learning engine |
| <p><b>Conclusions and suggestions</b><br/>AI-supported interpretation based on lung function: COPD</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 4. <b>Decision support:</b> most likely disease                                                    |
| <p>In addition to patient history and clinical exam, the following tests may help determine a final diagnosis:<br/>Repeat spirometry with bronchodilator test. Perform HRCT of the thorax.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 5. <b>Further suggestions:</b> best medical practice                                               |
| <p><b>Warnings</b><br/>The output of the report should always be considered in combination with patient history and clinical examination. The output is intended to support, not replace, clinical decision-making.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 6. <b>Warnings</b>                                                                                 |
| <p><b>Legend</b><br/>COPD Chronic Obstructive Pulmonary Disease<br/>OBD Other Obstructive Diseases (including: cystic fibrosis, bronchiectasis, bronchiolitis)<br/>Normal Normal lung function<br/>ILD Interstitial lung disease (including: idiopathic pulmonary fibrosis, nonspecific interstitial pneumonitis, sarcoidosis)<br/>NMD Neuromuscular disease (including: paralysis of the diaphragm, polymyositis, myopathy)<br/>PVD Pulmonary vascular disease (including: pulmonary hypertension, embolism, vasculitis)<br/>TD Thoracic deformity / Pleural disease (including: pneumectomy, lobectomy, chest wall problems, kyphoscoliosis)</p> <p>Physiological interpretation of the lung function tests: ERS/ATS 2021 (Stanojevic).<br/>Reference equations: Spirometry (GLI 2012), Static lung volumes (ECSS), Diffusion (GLI 2017 + 2020 correction).</p> <p>This report is approved for clinical use in the EU.<br/>Automatically generated by ArtiQ.PFT TEST-1fdca99 • Manufactured by ArtiQ NV • Leuven, Belgium</p> <p>info@ArtiQ.eu<br/>www.ArtiQ.eu</p> |                                                                                                    |

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QUNIQUE GmbH, Bahnhofweg 17, 5610 Wohlen, Schweiz