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
<p><b>ARTIQ</b></p> <p>Analyzed: 2025-10-27 13:11 Report ID: 699f20cd-9c3f-4508-a34a-9750adff5d60 Age: 64      Sex at birth: Female      Current Smoker: Yes      Pack-Years: 34</p>	
<p><b>Quality of lung function tests</b> FEV1 quality grade is A. FVC quality grade is B. The test quality is good and results can be interpreted with confidence.</p>	<p><b>1. Quality:</b> Quality of spirometry and meaning</p>
<p><b>Physiological interpretation of lung function tests</b> Moderate obstructive lung function. Bronchodilator response test is not performed. Hyperinflation and/or air trapping. Increased airway resistance. Moderate reduction of diffusion capacity.</p>	<p><b>2. Physiological interpretation:</b> spirometry description as per international standards</p>
<p><b>AI-supported interpretation</b></p>  <p>Legend for chart:  <span style="display:inline-block; width:10px; height:10px; background-color:#e0e0e0; border:1px solid #000;"></span> Asthma  <span style="display:inline-block; width:10px; height:10px; background-color:#c0c0c0; border:1px solid #000;"></span> COPD  <span style="display:inline-block; width:10px; height:10px; background-color:#a0a0a0; border:1px solid #000;"></span> OBD  <span style="display:inline-block; width:10px; height:10px; background-color:#808080; border:1px solid #000;"></span> Normal  <span style="display:inline-block; width:10px; height:10px; background-color:#606060; border:1px solid #000;"></span> ILD  <span style="display:inline-block; width:10px; height:10px; background-color:#404040; border:1px solid #000;"></span> NMD  <span style="display:inline-block; width:10px; height:10px; background-color:#202020; border:1px solid #000;"></span> PVD  <span style="display:inline-block; width:10px; height:10px; background-color:#000000; border:1px solid #000;"></span> TD</p>	<p><b>3. AI-supported interpretation:</b> disease probability estimation with machine learning engine</p>
<p><b>Conclusions and suggestions</b> AI-supported interpretation based on lung function: COPD</p>	<p><b>4. Decision support:</b> most likely disease</p>
<p>In addition to patient history and clinical exam, the following tests may help determine a final diagnosis: Repeat spirometry with bronchodilator test. Perform HRCT of the thorax.</p>	<p><b>5. Further suggestions:</b> best medical practice</p>
<p><b>Warnings</b> The output of the report should always be considered in combination with patient history and clinical examination. The output is intended to support, not replace, clinical decision-making.</p>	<p><b>6. Warnings</b></p>
<p><b>Legend</b></p> <p>COPD Chronic Obstructive Pulmonary Disease  OBD Other Obstructive Diseases (including: cystic fibrosis, bronchiectasis, bronchiolitis)  Normal Normal lung function  ILD Interstitial lung disease (including: idiopathic pulmonary fibrosis, nonspecific interstitial pneumonitis, sarcoidosis)  NMD Neuromuscular disease (including: paralysis of the diaphragm, poliomyelitis, myopathy)  PVD Pulmonary vascular disease (including: pulmonary hypertension, embolism, vasculitis)  TD Thoracic deformity / Pleural disease (including: pneumectomy, lobectomy, chest wall problems, kyphoscoliosis)</p> <p>Physiological interpretation of the lung function tests: ERS/ATS 2021 (Stanojevic).  Reference equations: Spirometry (GLI 2012), Static lung volumes (ECCS), Diffusion (GLI 2017 + 2020 correction).</p> <p>This report is approved for clinical use in the EU.  Automatically generated by ArtiQ.PFT TEST-1fdca99 • Manufactured by ArtiQ NV • Leuven, Belgium</p> <p style="text-align: right;">info@ArtiQ.eu www.ArtiQ.eu</p>	

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